mikrobiologie labor technik

Columbia Agar + 5% Sânge de Berbec

PRINCIPIU

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
Columbia Agar + 5% Sheep blood

PRINCIPLE

It is a non-selective medium for isolation and cultivation of microorganisms from clinical samples.

FORMULA

Components	g / l
Peptospecial	23.0
Starch	1.0
Sodium chloride	5.0
Sheep blood defibrinat	50.0 ml
Agar	14.0
Final pH 7.3 ± 0.2	

DESCRIPTION

Columbia Agar (with 5% sheep blood) is a non-selective medium containing blood for isolation and cultivation of demanding and tough to determine hemolytic reactions.

PROPERTIES

Peptospecial is a combination of peptones and is the source of nitrogen, carbon, sulfur and other factors essential for growth. Starch is used as an energy source. Sodium chloride maintains osmotic equilibrium of the environment. Sheep blood defibrinat provides additional factors of growth for organizations demanding and allow for the hemolytic reactions.

TECHNIQUE

Inoculate on the surface of plates with the sample loop using a sterile environment. Incubate the plates at $36 \pm 1^\circ\text{C}$ for 18-48 hours in anaerobic conditions or aerobically in the laboratory procedure.

INTERPRETATION OF RESULTS

Increases in microorganisms and hemolytic reactions.

There are 4 types of hemolysis:

- **α -haemolysis**: is the reduction of hemoglobin to methaemoglobin in the environment around the colony, and causing green discoloration of the environment;
- **β -haemolysis**: is lysis of erythrocytes, which show a clear zone around the colony;
- **γ -haemolysis**: no red blood cells are destroyed and do not notice any change in the culture medium;
- **α' -haemolysis**: has a small area of complete hemolysis surrounded by a zone of partial lysis.

BIBLIOGRAPHY.

1. Ellner, P.D., C.J. Stoessel., E. Drakeford, and F. Vasi (1966). A new culture medium for medical bacteriology. *Am. J.Clin. Path.* 45, 502-504.
2. Isenberg, H.D. (ed.) (1992). *Clinical microbiology procedures handbook*, vol. 1 American Society for Microbiology, Washington, D

PACKAGING

Dehydrated medium

To be stored between 1 and 30°C

140121A: Flask of 500 g

140121L: Flask of 500 g

Ready poured medium + Sheep's blood

To be stored between 2 and 8°C

110122: Pack of 20 dishes 90 mm \varnothing

120179: Pack of 32 tubes of 7 ml