mikrobiologie labor technik

A.D.C.L. (HYNES)

PRINCIPIU

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
A.D.C.L. (HYNES)

PRINCIPLE

A.D.C.L. (Desoxycholate-Citrate-Lactose-Agar) is a selective medium for isolation of *Salmonella* and *Shigella*. High levels of desoxycholate and citrate allow to inhibit the growth of Gram positive bacteria and partially those of a lot of Enterobacteria. The fermentation of lactose induces an acidification and a colour change of the neutral red. The reduction of thiosulfate is revealed by ferric ions in black iron sulphide form.

FORMULA

Component	g/l
Polypeptone	5.00
Meat extract	5.00
Sodium citrate	8.50
Sodium desoxycholate	5.00
Lactose	10.00
Ferric ammonium citrate	1.00
Neutral red	0.02
Sodium thiosulfate	5.40
Agar	12.00
Final pH : 7.4 ± 0.2 at 25°C	

PREPARATION

Pour 52 g of powder in one litre of purified water. Bring slowly to the boil under agitation until complete dissolution. Maintain the boiling during one minute. **DO NOT AUTOCLAVE.**

PROCEDURE

When the medium is liquefied, cool to 47 ± 2°C and dispatch into sterile Petri plates.

Culture of stools :

Inoculate 2 plates from suspensions prepared with stools. Inoculate at the same time sodium selenite broth and Muller Kauffmann medium and incubate for 24 hours at 37°C.

Detection of *Salmonella* in pharmaceutical products

Inoculate 10 g of product to be tested in 100 ml of lactose broth for 5 to 24 hours at 37°C. Take 10 ml of this pre-culture and inoculate 100 ml of T.B.G. broth. Incubate 18-24 hours at 42-43°C and realize subcultures on 2 media among : D.C.L., X.L.D. and V.B.R.P. Kristensen.

RESULTS

Note the appearance of the colonies.

Escherichia coli : pink to red, 1 to 2 mm diameter, with halo of precipitation (nevertheless, most of strains are inhibited).

Shigella sonnei : colourless to pale pink (low degradation of lactose).

Shigella flexneri : colourless to pale pink, with rounded centre.

Salmonella Paratyphi B : colourless to pale pink after 18 hours, slightly opaque with black centre after 48 hours.

Salmonella Typhi : colourless to pale pink after 18

hours, slightly opaque with grey centre after 48 hours.

Others *Salmonella* H₂S + : equal to *S. Paratyphi B*.

***Proteus* :** equal to *Salmonella* or *Shigella*, brilliant with black centre and characteristic odour.

LIMITS AND PRECAUTIONS

Lactose + micro-organisms grow as pink colonies with a halo of precipitation of desoxycholic acid.

Colonies of lactose - micro-organisms are colourless. This medium may precipitate

with the cold. So it is advised not to store it between 2 and 8°C.

BIBLIOGRAPHY

1. Hynes M. 1942. The isolation of intestinal pathogens by selective media. J. Pathol. Bacteriol. 54:193-207.

PACKAGING

Dehydrated medium

(To be stored between 1 and 30°C)

140101A: Flask of 500 g

Ready to use medium

(To be stored between 15 and 25°C)

130101: Pack of 3 flasks of 200 ml

110101: Pack of 20 plates of 90 mm Ø

120175: Pack of 32 tubes of 7 ml