mikrobiologie labor technik

Mac Conkey Agar Nr. 3

PRINCIPIU

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2. Mac Conkey A. 1908. Bile salt media and their advantage in some bacteriological examination. J. Hyg.,8:322-334.
3. Pharmacopée Européenne. Milieu H.

AMBALARE

Mediul deshidratat

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140139A: Flacon de 500 g

140139L: Flacon de 500 g


Mediul gata preparat

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110139: Cutie cu 20 de plăci de 90 mm Ø

130139: Cutie cu 3 sticle de 200 ml

120183: Cutie de 32 tuburi de 7 ml

mikrobiologie labor technik

Mac Conkey Agar N°3

PRINCIPLE

Mac Conkey agar (Mac Conkey N°3) is a selective and differential plating medium mainly used for the detection and isolation of gram-negative organism from dairy, foodstuff and waters samples. In clinical specimen this medium is adapted for the screening of pathogen enterobacteria especially in babies stools. Pharmaceutical and cosmetic industries this medium is used to detect the presence or the absence of *Escherichia coli*. Bile salts and Crystal Violet inhibits Gram positive germs growth. The pH indicator, neutral red, helps to differentiate lactose-fermenting from lactose none fermenting gram negative enteric bacilli.

FORMULA

Components	g/l
Peptone	17.00
Proteose Peptone	3.00
Monohydrated Lactose	10.00
Bile salts	1.50
Sodium chloride	5.00
Neutral red	0.03
Crystal violet	0.001
Agar	13.50
Final pH : 7.1 ± 0.2 at 25°C	

METHOD

Pour 50 g of powder into 1 litre purified water. Bring slowly to the boil, homogenise until complete dissolution. Sterilize in autoclave for 15 minutes at 121°C. Liquefy the medium then cool to 45-50°C and dispense into sterile Petri plates, tubes or flasks. Inoculate the surface in order to obtain isolated colonies.

PROCEDURE

Application n°1:

When analysing, water, dairy, foodstuff, dispense 1 ml of the sample or its decimal dilutions into a sterile Petri plate then pour about 15 ml of medium (45-57°C), homogenise and let the medium set. Incubate 18 to 24 hours at 37°C

Application n°2:

When screening for pathogen *Enterobacteria* in babies stools or urine swab the specimen (or the enriched specimen) onto the surface of a Mac Conkey plate. Repeat the procedure with a second selective medium such as D.C.L.S for *Salmonella* and *Shigella*. Incubate 18 to 24 hours at 37°C.

Application n°3:

When screening for *Escherichia coli* in pharmaceutical and cosmetic products subculture the enriched sample Mac Conkey broth for pharmaceutical preparations and LT100 for cosmetic products) onto the surface of a Mac Conkey plate. Incubate at 32.5 ± 2.5°C for 18 to 72 hours.

RESULTS

Lactose – germs grow as colourless colonies as oppose to lactose + germs that grow as pink to brick- red colonies with or with out a precipitation zone.

BIBLIOGRAPHY

1. Mac Conkey A. 1905. Lactose-fermenting bacteria in faeces. J. Hyg. 8:333-379.
2. Mac Conkey A. 1908. Bile salt media and their advantage in some

bacteriological examination. J. Hyg.,
8:322-334.

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PACKAGING

Dehydrated medium (Store between 1 and 30°C)

140139A: Flask of 500 g

140139L: Flask of 500 g

Ready to use medium(Store between 2 and 25°C)

130139: Pack of 3 flasks of 200 ml

110139: Pack of 20 plates 90 mm Ø

120183: Pack of 32 tubes of 7 ml