mikrobiologie labor technik

Selenit-Cistină Bulion

PRINCIPIU

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
Selenit-Cystine Broth

PRINCIPLE

Selenite Cystine Broth is a selective enrichment broth used for the detection of *Salmonella* in stools, foodstuffs and pharmaceutical products. Sodium Selenite helps to slow and/or inhibit the growth of coliform and intestinal streptococci.

FORMULA (single concentrated broth)

Components	g/l
Tryptone	5.00
Lactose	4.00
Disodium phosphate (12 H ₂ O)	10.00
Sodium Selenite	4.00
L-Cystine	0.01
Final pH : 7.0 ± 0.2 at 25°C	

PREPARATION

Suspend 23 g of powder per litre of purified water. Bring slowly to the boil as to completely dissolve the powder boil for one minute. **DO NOT AUTOCLAVE**. Dispatch in tubes or flask

METHOD

Cool the medium to 25°C. Inoculate 10 ml of sample into 100 ml of broth (or 1 ml in 10 ml). According to the origin of the sample, incubate at 37 or 43°C for 24 hours.

RESULTS

After incubation, subculture onto different selective Agar for *Salmonella*.

LIMITS AND PRECAUTIONS

For detection of specified germs in urine samples, this broth should be prepared double concentrated and inoculated with the equal parts of urine and broth. Over incubating the tubes will not secure a good recovery of *S. typhi* & *paratyphi* since these germs develop faster than coliforms. Prepared medium must be kept at 2-8°C in the dark and used within a week. One should not use tubes or flasks that show the presence of a red precipitate in the bottom. On the contrary an orange precipitate will not harm the performances of the broth.

BIBLIOGRAPHY

1. FIL93. 1980. Norme Internationale provisoire.
2. Leifson E. 1936. New Selenite Enrichment Media for the isolation of Typhoid and Paratyphoid (*Salmonella*) Bacilli. Am. J. Hyg. 24:423-432.
3. USP XXIII (1995).

PACKAGING

Dehydrated (Store between 1 and 30°C)
140198A: Flask of 500 g
140198L: Flask of 500 g

Ready to use (Store between 2 and 8°C)
Single concentration
130225: Pack of 3 flasks of 200 ml
120225: Pack of 32 tubes of 7 ml